



*Hygienic & Clean
Air Sampler*




HYDRA SAMPLER

HYgienic Distributed Rotation Air sampler


► **Up to 8 Air Samplers can be controlled**

**Continuous
sampling**



*Validated according to EN 17141
and ISO 14698 by an external
independent laboratory*

► Continuous monitoring with 3 color codes





Ready to sample




Remove plate




Sampling in progress




1. Insert plates in all air sampler



2. Select and run the first sampler




3. At sampling completion, start second sampler



4. Proceed with the rotation, reset the sampler when plate is replaced

► Bidirectional communication



- Independent sampling units
- Automatic identification of sampler serial number
- Integrated air flow sensor, the length of the cable does not affect the airflow

► Data integrity




Data can be exported to:

- Secured PDF file onto USB stick
- External printer
- DCS network

21 CFR PART 11 version available

► High quality stainless steel connectors





Via S. Predengo, 27/29
26022 - Castelverde | Cremona | ITALY

www.hyclas.com